

HAND ACTIVATED MONOPOLAR FORCEPS



Fine tips for accuracy and speed of tissue dissection



Hand-activated switch eliminates need for foot pedals



Attached cord can be used with most ESUs

- Fine point, accurate, electrosurgical tissue forceps
- New standard for tissue dissection and hemostasis
- Available in two lengths for body and face procedures

B 65245	7"	17.5 cm	Serrated Jaw
B 65246	7"	17.5 cm	Smooth Jaw
B 65247	8"	20 cm	Smooth Jaw, Micro-Tip
B 65244	10"	25.5 cm	Serrated Jaw
B 65242	10"	25.5 cm	Smooth Jaw



Serrated Jaw



Smooth Jaw



Smooth Jaw, Micro-Tip