


FINE INSTRUMENTS AND EQUIPMENT
FOR LIPOSUCTION AND FAT GRAFTING



LIPOSUCTION CANNULAS

HANDLES AND HUBS

Vitruvian™ Handle



Luer Lock Hub



Luer Lock Aspiration Handle



B 89ASH-LL

60cc Toomey Syringe Hub



Toomey Aspiration "Delrin" Handle



B 89TAH-60CC

Power Aspiration Hub



CANNULA PATTERNS



Mercedes



Double Mercedes




Accelerator



Mirrored Accelerator



12-Hole Large Volume Harvester



Candy Cane



Helix



Sattler



Spatula



Basket



Double Basket



Fat Disruptor



Standard

POPULAR PATTERNS

Accelerator

B 8III-215L	2mm X 15cm
B 8III-215-SYR	2mm X 15cm
B 8III-226L	2mm X 26cm
B 8III-315L	3mm X 15cm
B 8III-326L	3mm X 26cm
B 8III-326-SYR	3mm X 26cm
B 8III-332L	3mm X 32cm
B 8III-336L	3mm X 36cm
B 8III-315-SYR	3mm X 15cm
B 8III-415L	4mm X 15cm
B 8III-415-SYR	4mm X 15cm
B 8III-426L	4mm X 26cm
B 8III-426-SYR	4mm X 26cm
B 8III-432L	4mm X 32cm
B 8III-432-SYR	4mm X 32cm
B 8III-436L	4mm X 36cm
B 8III-526L	5mm X 26cm
B 8III-532L	5mm X 32cm

Basket

B 8BAS-326L	3mm X 26cm
B 8BAS-426L	4mm X 26cm
B 8BAS-432L	4mm X 32cm
B 8BAS-526L	5mm X 26cm
B 8BAS-532L	5mm X 32cm

Double Basket

B 8DBAS-426L	4mm X 26cm
B 8DBAS-432-SYR	4mm X 32cm
B 8DBAS-526L	5mm X 26cm

Candy Cane

B 8CC-326L	3mm X 26cm
B 8CC-332L	3mm X 32cm
B 8CC-426L	4mm X 26cm
B 8CC-432L	4mm X 32cm
B 8CC-436L	4mm X 36cm
B 8CC-526L	5mm X 26cm
B 8CC-532L	5mm X 32cm

Fat Disruptor

B 8MBAT-315L	3mm X 15cm
B 8MBAT-326L	3mm X 26cm
B 8MBAT-426L	4mm X 26cm
B 8MBAT-432L	4mm X 32cm
B 8MBAT-532L	5mm X 32cm

Large Volume Harvester

B 8LVH-426L	4mm X 26cm
B 8LVH-432L	4mm X 32cm

Mercedes

B 8MER-215L	2mm X 15cm
B 8MER-226L	2mm X 26cm
B 8MER-315L	3mm X 15cm
B 8MER-315-SYR	3mm X 15cm
B 8MER-326L	3mm X 26cm
B 8MER-326-SYR	3mm X 26cm
B 8MER-332L	3mm X 32cm
B 8MER-336L	3mm X 36cm
B 8MER-415L	4mm X 15cm
B 8MER-415-SYR	4mm X 15cm
B 8MER-426L	4mm X 26cm
B 8MER-426-SYR	4mm X 26cm
B 8MER-432L	4mm X 32cm
B 8MER-432-SYR	4mm X 32cm
B 8MER-436L	4mm X 36cm
B 8MER-526L	5mm X 26cm
B 8MER-532L	5mm X 32cm
B 8MER-536L	5mm X 36cm

Sattler

B 8SAT-215L	2mm X 15cm
B 8SAT-226L	2mm X 26cm
B 8SAT-315L	3mm X 15cm
B 8SAT-326L	3mm X 26cm
B 8SAT-332L	3mm X 32cm

Spatula

B 8SPA-215L	2mm X 15cm
B 8SPA-315L	3mm X 15cm
B 8SPA-415L	4mm X 15cm
B 8SPA-515L	5mm X 15cm
B 8SPA-615L	6mm X 15cm

Standard 1

B 8STD-315L	3mm X 15cm
B 8STD-315-SYR	3mm X 15cm
B 8STD-426L	4mm X 26cm
B 8STD-526L	5mm X 26cm

CANNULA ORDERING INFORMATION

Please use the following steps in making your selection:

- Choose a pattern (i.e. Mercedes)
- Choose a diameter (i.e. 3mm)
- Choose a length (i.e. 26cm)
- Choose a handle design (i.e. large handle)


Gauge (ga):	20	18	16	14	12	11	8	6
Diameter (mm):	0.9	1.2	1.6	2.0	2.7	3.0	4.0	5.0

POWER LIPOSUCTION CANNULAS

MICROAIRE®-TYPE POWER LIPOSUCTION CANNULAS

Using tight tolerances to minimize vibration, this superior, reusable cannula is fabricated from the highest-quality materials.

- Reusable / Autoclavable
- Cannula tubing connector designed for a non-slip fit for both aspiration and infiltration
- Full line of sizes and patterns
- High quality materials and craftsmanship
- Made in the USA from patented design
US Patent D950,045 / D950,046



POWER LIPOSUCTION STARTER SET

B 8B-MER-326 PWR	Mercedes Power Cannula 3mm × 26cm
B 8B-MER-432 PWR	Mercedes Power Cannula 4mm × 32cm
B 8B-BAS-432 PWR	Basket Power Cannula 4mm × 32cm
B 8B-III-326 PWR	Accelerator Power Cannula 3mm × 26cm
B 8B-III-432 PWR	Accelerator Power Cannula 4mm × 32cm

POPULAR PATTERNS - POWER LIPOSUCTION CANNULAS



ACCELERATOR

8B-III-315 PWR	3MM × 15CM STRAIGHT
8B-III-326 PWR	3MM × 26CM STRAIGHT
8B-III-332 PWR	3MM × 32CM STRAIGHT
8B-III-426 PWR	4MM × 26CM STRAIGHT
8B-III-432 PWR	4MM × 32CM STRAIGHT
8B-III-436 PWR	4MM × 36CM STRAIGHT
8B-III-526 PWR	5MM × 26CM STRAIGHT
8B-III-532 PWR	5MM × 32CM STRAIGHT



MIRRORED ACCELERATOR

8B-MIII-326 PWR	3MM × 26CM STRAIGHT
8B-MIII-332 PWR	3MM × 32CM STRAIGHT
8B-MIII-426 PWR	4MM × 26CM STRAIGHT
8B-MIII-432 PWR	4MM × 32CM STRAIGHT
8B-MIII-526 PWR	5MM × 26CM STRAIGHT
8B-MIII-532 PWR	5MM × 32CM STRAIGHT



BASKET

8B-BAS-326 PWR	3MM × 26CM STRAIGHT
8B-BAS-332 PWR	3MM × 32CM STRAIGHT
8B-BAS-426 PWR	4MM × 26CM STRAIGHT
8B-BAS-426 ANG PWR	4MM × 26CM ANGLED
8B-BAS-432 PWR	4MM × 32CM STRAIGHT
8B-BAS-432 ANG PWR	4MM × 32CM ANGLED
8B-BAS-440 PWR	4MM × 40CM STRAIGHT
8B-BAS-440 ANG PWR	4MM × 40CM ANGLED
8B-BAS-526 PWR	5MM × 26CM STRAIGHT
8B-BAS-526 ANG PWR	5MM × 26CM ANGLED
8B-BAS-532 PWR	5MM × 32CM STRAIGHT
8B-BAS-532 ANG PWR	5MM × 32CM ANGLED
8B-BAS-536 PWR	5MM × 36CM STRAIGHT
8B-BAS-536 ANG PWR	5MM × 36CM ANGLED



DOUBLE BASKET

8B-DBAS-326 PWR	3MM × 26CM STRAIGHT
8B-DBAS-332 PWR	3MM × 32CM STRAIGHT
8B-DBAS-332 ANG PWR	3MM × 32CM ANGLED
8B-DBAS-426 PWR	4MM × 26CM STRAIGHT
8B-DBAS-432 PWR	4MM × 32CM STRAIGHT
8B-DBAS-432 ANG PWR	4MM × 32CM ANGLED
8B-DBAS-440 PWR	4MM × 40CM STRAIGHT
8B-DBAS-440 ANG PWR	4MM × 40CM ANGLED
8B-DBAS-526 PWR	5MM × 26CM STRAIGHT
8B-DBAS-532 PWR	5MM × 32CM STRAIGHT
8B-DBAS-532 ANG PWR	5MM × 32CM ANGLED
8B-DBAS-540 PWR	5MM × 40CM STRAIGHT
8B-DBAS-540 ANG PWR	5MM × 40CM ANGLED



SPATULA

8B-SPA-315 PWR	3MM × 15CM STRAIGHT
8B-SPA-415 PWR	4MM × 15CM STRAIGHT



MERCEDES

8B-MER-215 PWR	2MM × 15CM STRAIGHT
8B-MER-315 PWR	3MM × 15CM STRAIGHT
8B-MER-326 PWR	3MM × 26CM STRAIGHT
8B-MER-332 PWR	3MM × 32CM STRAIGHT
8B-MER-332 ANG PWR	3MM × 32CM ANGLED
8B-MER-426 PWR	4MM × 26CM STRAIGHT
8B-MER-426 ANG PWR	4MM × 26CM ANGLED
8B-MER-432 PWR	4MM × 32CM STRAIGHT
8B-MER-432 ANG PWR	4MM × 32CM ANGLED
8B-MER-436 PWR	4MM × 36CM STRAIGHT
8B-MER-436 ANG PWR	4MM × 36CM ANGLED
8B-MER-526 PWR	5MM × 26CM STRAIGHT
8B-MER-526 ANG PWR	5MM × 26CM ANGLED
8B-MER-532 PWR	5MM × 32CM STRAIGHT
8B-MER-532 ANG PWR	5MM × 32CM ANGLED
8B-MER-536 PWR	5MM × 36CM STRAIGHT
8B-MER-536 ANG PWR	5MM × 36CM ANGLED



DOUBLE MERCEDES

8B-DMER-326 PWR	3MM × 26CM STRAIGHT
8B-DMER-426 PWR	4MM × 26CM STRAIGHT
8B-DMER-432 PWR	4MM × 32CM STRAIGHT
8B-DMER-436 PWR	4MM × 36CM STRAIGHT
8B-DMER-526 PWR	5MM × 26CM STRAIGHT
8B-DMER-532 PWR	5MM × 32CM STRAIGHT
8B-DMER-536 PWR	5MM × 36CM STRAIGHT



HELIX 3-PORT

8B-HLX-326 PWR	3MM × 26CM STRAIGHT
8B-HLX-332 PWR	3MM × 32CM STRAIGHT
8B-HLX-426 PWR	4MM × 26CM STRAIGHT
8B-HLX-432 PWR	4MM × 32CM STRAIGHT
8B-HLX-436 PWR	4MM × 36CM STRAIGHT
8B-HLX-526 PWR	5MM × 26CM STRAIGHT
8B-HLX-532 PWR	5MM × 32CM STRAIGHT
8B-HLX-536 PWR	5MM × 36CM STRAIGHT



LARGE VOLUME HARVESTER

8B-LVH-326 PWR	3MM × 26CM STRAIGHT
8B-LVH-332 PWR	3MM × 32CM STRAIGHT
8B-LVH-426 PWR	4MM × 26CM STRAIGHT
8B-LVH-432 PWR	4MM × 32CM STRAIGHT
8B-LVH-440 PWR	4MM × 40CM STRAIGHT
8B-LVH-526 PWR	5MM × 26CM STRAIGHT
8B-LVH-532 PWR	5MM × 32CM STRAIGHT



CANDY CANE

8B-CC-326 PWR	3MM × 26CM STRAIGHT
8B-CC-426 PWR	4MM × 26CM STRAIGHT
8B-CC-432 PWR	4MM × 32CM STRAIGHT
8B-CC-436 PWR	4MM × 36CM STRAIGHT
8B-CC-526 PWR	5MM × 26CM STRAIGHT
8B-CC-532 PWR	5MM × 32CM STRAIGHT

TUMESCENT INFILTRATION

B 89055

VITRUVIAN™ VARIABLE SPEED HIGH FLOW INFILTRATION PUMP WITH VOLUME READ OUT

The Black & Black Vitruvian™ Infiltration Pump is a counting peristaltic pump that uses three rollers to create a smooth pumping action to deliver tumescent fluid to the operative site.

Powered by a state-of-the-art stepper motor, this machine provides a smooth and consistent pumping action at a rate of up to 1000mL/min.

The Vitruvian™ Infiltration Pump was designed with industry leading features such as fluid counting technology, a rate adjusting knob, prime button, reverse switch, and dual foot pedals that can be used in “momentary” or “maintained” mode. Large, bright LEDs are easily viewed from anywhere in the room and clearly display both volume and rate.

- Peristaltic pump with three rollers
- Dual-voltage motor – 100V/240V
- Bright, easy-to-read display screen
- Volume [mL] and Rate [mL/min] displayed
- Reverse switch for flow direction preference
- Prime button for fast filling of tubing line
- Momentary / maintained pneumatic foot-switch functions (For ON/OFF action or continuous action, respectively)
- Dual foot pedals



INFILTRATION TUBING *For use with the Vitruvian™ Infiltration Pump*



B 89213

Large Bore™ SofTouch™ Pump Tubing, single spike, 13' length, 10 per box. Single use, sterile.



B 89210

Large Bore™ SofTouch™ Tubing with molded connectors for pump systems. 13' length, 10 per box. Single use, sterile.



B 89214

Dual Spike Large Bore™ SofTouch™ tubing for infiltration pump systems. 13' length, 10 per box. Single use, sterile.

TUMESCENT INFILTRATION CANNULAS

Body



B 8TI-1015-LL	10ga X 15cm
B 8TI-1032-LL	10ga X 32cm
B 8TI-1215-LL	12ga X 15cm
B 8TI-1220-LL	12ga X 20cm
B 8TI-1226-LL	12ga X 26cm
B 8TI-1232-LL	12ga X 32cm
B 8TI-1236-LL	12ga X 36cm
B 8TI-1415-LL	14ga X 15cm
B 8TI-1420-LL	14ga X 20cm
B 8TI-1426-LL	14ga X 26cm
B 8TI-1432-LL	14ga X 32cm
B 8TI-1436-LL	14ga X 36cm
B 8TI-1615-LL	16ga X 15cm
B 8TI-1620-LL	16ga X 20cm
B 8TI-1626-LL	16ga X 26cm
B 8TI-1632-LL	16ga X 32cm

Facial



B 8FAN-1415-LL	14ga X 15cm
B 8FAN-1615-LL	16ga X 15cm

LL = Luer Lock

HANDLES



B 89TIH-C
Infiltration Handle With On/Off Control



B 89IH-NC
Infiltration Handle Without Controls



B 89IH-C
Infiltration Handle With On/Off Control

INFILTRATION TUBING

For use with non-Vitruvian™ Infiltration Pumps



B 89216

Large Bore™ Tubing with silicone insert for pump systems.
13' length, 10 per box.
Single use, sterile.



B 89215

Dual Spike Infiltration Tubing with silicone insert for pump systems.
16' length, 10 per box.
Single use, sterile.



B 89217

Large Bore™ Transfer Tubing with silicone insert for pump systems.
Single spike, 16' length, 10 per box.
Single use, sterile.

ASPIRATOR SYSTEM

B 89025

VITRUVIAN™ ULTIMATE ASPIRATOR

The Black & Black Vitruvian™ Ultimate Aspirator is one of the most powerful, quiet and versatile liposuction aspirators on the market. Powered by two oil-less, rocking piston vacuum pumps, this machine provides excellent suction power along with quiet operation. Its dual voltage motor can be used with 120V/230V electrical service.

Our Vitruvian™ Ultimate Aspirator is carefully designed for optimal performance, with an upright stance and shopping-cart style handle for ease of mobility. Bright LED vacuum gauges are easy to read. The collection canisters are mounted for ease of viewing during surgery. The single knob control can be adjusted 0 to 29 in. of Hg vacuum.

- Dual pumps for optimal performance – 1 HP total
- Flow rate – 6.4 cubic ft./min
- Quiet operation (50 dB)
- Dual pneumatic foot-switches
- Large 3" hospital-grade wheels
- Bright, easy-to-read vacuum gauges
- Easy-to-see collection canisters
- Extendable pole to hold fluid bags
- Dual-voltage motors – 120V/230V



HIGH FLOW ASPIRATION TUBING



B 89200

SoftTouch™ formulation provides excellent drapability.
10' length, 10 per box.
Single use, sterile.



B 89205

Tapered Aspiration Tubing for power-assisted liposuction. SoftTouch™ formulation provides excellent drapability.
10' length, 10 per box.
Single use, sterile.

COLLECTION CANISTERS AND LINERS



Vitruvian™ Collection Canister

Designed for high vacuum liposuction procedures and made with lipid-resistant materials.

B 89140

- 2,000 CC / mL Canisters – 12 per case
- Disposable, non-sterile
- Excellent clarity for enhanced visualization of contents
- Lids contain multiple ports for versatility and tandem set up
- Integrated Aerostat® filter with >99.99% bacterial efficiency
- Hi-Flow shut off valve closes automatically to prevent overflow and vacuum line contamination
- Vacuum level rating of 29" Hg/736mm Hg
- Nestable for storage convenience
- Latex, BPA, and DEHP free



Hospira (Abbott, IP) Reusable Canister

B 89124

- 2,000 ml

B 89125

Disposable Liners

- 2,000 ml
- 10 per box



Baxter Reusable Canister

B 89130

- 1,500 ml

B 89135

Disposable Liners

- 1,500 ml
- 10 per box



B 89104

Aspirator filter 0.3 micron for Vitruvian™ Ultimate Aspirator

B 89105

0.3 micron filter for most other aspirators



B 89220

Filter Connecting Tube

- 30" non-sterile
- 10 per box



B 89026A

Pneumatic Foot Pedal for Aspirator

B 89026P

Pneumatic Foot Pedal for Infiltration Pump



B 89030

20AMP Hospital Grade Power Cord 10ft length for Aspirator

B 89032

15AMP Hospital Grade Power Cord 10ft length for Pump

B 89031


20AMP Hospital Grade Power Strip 25ft length

FAT GRAFTING

TISSUE INJECTION CANNULAS



Curve Options



Micro Tissue Injection

Type I

B 8FIC-I-167-LL	16ga X 7cm
B 8FIC-I-169-LL	16ga X 9cm
B 8FIC-I-187-LL	18ga X 7cm
B 8FIC-I-189-LL	18ga X 9cm
B 8FIC-I-195-LL	19ga X 5cm w/ Stylette
B 8FIC-I-203-LL	20ga X 3cm w/ Stylette
B 8FIC-I-205-LL	20ga X 5cm w/ Stylette
B 8FIC-I-CC1610-LL	16ga X 10cm, Concave
B 8FIC-I-CV1610-LL	16ga X 10cm, Convex

Type II

B 8FIC-II-167-LL	16ga X 7cm
B 8FIC-II-169-LL	16ga X 9cm
B 8FIC-II-187-LL	18ga X 7cm
B 8FIC-II-189-LL	18ga X 9cm
B 8FIC-II-CC1610-LL	16ga X 10cm, Concave

Type III

B 8FIC-III-167-LL	16ga X 7cm
B 8FIC-III-169-LL	16ga X 9cm
B 8FIC-III-187-LL	18ga X 7cm
B 8FIC-III-189-LL	18ga X 9cm

Body Tissue Injection

Type I

B 8FIC-I-1215-LL	12ga X 15cm
B 8FIC-I-1220-LL	12ga X 20cm
B 8FIC-I-1225-LL	12ga X 25cm
B 8FIC-I-1415-LL	14ga X 15cm
B 8FIC-I-1420-LL	14ga X 20cm
B 8FIC-I-1425-LL	14ga X 25cm

V-Tip Tissue Injection

B 8FIC-V-167-LL	16ga X 7cm
B 8FIC-V-169-LL	16ga X 9cm
B 8FIC-V-187-LL	18ga X 7cm

Toledo V-Tip Tissue Injection

B 8FIC-V-1215-LL	12ga X 15cm
B 8FIC-V-1225-LL	12ga X 25cm
B 8FIC-V-1415-LL	14ga X 15cm
B 8FIC-V-149-LL	14ga X 9cm

Fat Transfer



Luer-to-luer connector

B 89LL2	2mm
B 89LL3	3mm
B 89LL4	4mm
B 89LL2.4	2.4mm
B 89LL1.4	1.4mm
B 89LL1.2	1.2mm



Becker luer-to-luer morselizer

B 89BEC-LL3

MICRO ASPIRATION CANNULAS

Accelerator



B 8III-1215-LL	12ga X 15cm
B 8III-1225-LL	12ga X 25cm
B 8III-1415-LL	14ga X 15cm

FAC



B 8FAC-1215-LL	12ga X 15cm
B 8FAC-1215-LL-CVD	12ga X 15cm
B 8FAC-1220-LL	12ga X 20cm
B 8FAC-1225-LL	12ga X 25cm
B 8FAC-1415-LL	14ga X 15cm
B 8FAC-1425-LL	14ga X 25cm

Fat Disruptor



B 8MBAT-1215-LL	12ga X 15cm
B 8MBAT-1225-LL	12ga X 25cm
B 8MBAT-1415-LL	14ga X 15cm

Large Volume Harvester



B 8LVH-1220-LL	12ga X 20cm
----------------	-------------

Mercedes



B 8MER-1215-LL	12ga X 15cm
B 8MER-1225-LL	12ga X 25cm
B 8MER-1415-LL	14ga X 15cm

Sattler



B 8SAT-1215-LL	12ga X 15cm
B 8SAT-1225-LL	12ga X 25cm
B 8SAT-1415-LL	14ga X 15cm
B 8SAT-1425-LL	14ga X 25cm

FAT COLLECTION CANISTER KITS

Canister Kits are completely autoclavable and include:

- Reusable canister, lid and stand
- Large bore universal taper fitting
- 6' silicone suction tubing with coupling piece
- Silicone transfer tube with C-Clamp + Luer fitting
- Silicone transfer tube with C-Clamp + Toomey syringe fitting
- Silicone Transfer Tube With C-Clamp + Tapered Port fitting

B 89110	1000ml Canister kit
B 89116	2000ml Canister kit
B 89117	3000ml Canister kit



RE-ORDER PARTS

B 89107	Silicone transfer tube w/Luer Lock fitting w/clamp
B 89108	Stand
B 89109	Silicone transfer tube w/Toomey syringe fitting w/clamp
B 89119	Lid
B 89120	Coupling piece for canister lid
B 89121	6' silicone suction tubing w/coupling piece

ACCESSORIES

SYRINGE LOCKS

B 89SL-10CC	Snap Lock For 10cc Syringe
B 89SL-20CC	Snap Lock For 20cc Syringe
B 89SL-30CC	Snap Lock For 30cc Syringe
B 89SL-60CC	Snap Lock For 60cc Syringe



Snap Locks

B 89TL-10CC	Twist Lock "Horseshoe" 10cc Syringe
B 89TL-20CC	Twist Lock "Horseshoe" 20cc Syringe
B 89TL-30CC	Twist Lock "Horseshoe" 30cc Syringe
B 89TL-60CC	Twist Lock "Horseshoe" 60cc Syringe



Twist Lock "Horseshoes"

SYRINGES

B 89SYR-60-TT	60cc Monoject Syringe With Toomey Tip, Sterile, 20 Per Box
B 89SYR-60-LL	60cc Luer Lock Syringe, Sterile, 40 Per Box
B 89SYR-10-LL	10cc Luer Lock Syringe, Sterile, 100 Per Box
B 89SYR-1-LL	1cc Luer Lock Syringe, Sterile, 100 Per Box




60cc Monoject Syringe With Toomey Tip

B 89PCG-1000	Pressure Cuff Infiltration System With Gauge For 1000ml Bag
--------------	---

SYRINGE RACKS

B 89SR-10

Syringe rack for 1cc – 10cc, autoclavable



B 89SR-60

Syringe rack for 60cc, autoclavable



CLEANING BRUSHES

- | | |
|-------------|------------------------------------|
| B 89CLB-116 | Cleaning Brush 1.6mm x 39cm (16ga) |
| B 89CLB-123 | Cleaning Brush 2mm x 60cm |
| B 89CLB-145 | Cleaning Brush 3mm - 4mm x 60cm |
| B 89CLB-156 | Cleaning Brush 5mm - 6mm x 60cm |



SYRINGE-TO-LUER



B 89TML
Toomey 60cc to Luer Lock
Cannula Connector
Stainless Steel



B 89TML
Toomey 60cc to Luer Lock
Syringe Connector
Stainless Steel



B 89605

Entrese Variable Depth Punch

ACCESSORIES

LipoPort™



LipoPort™ is a must-have for any surgeon performing liposuction. It improves the aesthetic outcome of incisions and aids in wound healing by saving the skin from friction burns.

LipoPort™ can be used with both manual and power-assisted liposuction procedures. This product is reusable and made in the USA.

B 89LP34


Medium: accommodates
2–4mm cannula

B 89LP56

Large: accommodates
2–5mm basket cannula and
up to 6mm standard cannula

- Saves the skin from friction burns and aids in wound healing
- Improves aesthetic outcome of incisions
- Reusable /autoclavable
- Easy to insert and anchor with 3-0 nylon suture
- Accommodates straight and angled power-assisted, laser-assisted, and standard cannulas
- External notch facilitates optimal cannula orientation
- Two convenient sizes that accommodate most cannulas

ContourFoam™



ContourFoam™, designed to mimic the original compression foam Reston™, allows for uniform compression after liposuction when used with compression garments.

- Uniform compression
- Reduces bruising
- Reduces swelling
- Large size roll option for seam-free application
- Flexible polyurethane foam for worry-free use
- Economical

Pre-Cut ContourFoam Sheets

B 89601 .5" x 14" x 9" sheet; 50/box



ContourFoam Rolls

B 89600 .5" x 14" x 72" roll; 4/box

B 89602 .5" x 14" x 48" roll; 4/box



POPULAR LIPOPLASTY SETS

AUTOLOGOUS FAT TRANSFER SET

B 8FAC-1215-LL	FAC 12ga X 15cm luer lock harvesting cannula
B 8FIC-I-167-LL	Type I 16ga X 7cm fat infiltration cannula
B 8FIC-II-169-LL	Type II 16ga X 9cm fat infiltration cannula
B 8SAT-1415-LL	Sattler 14ga X 15cm luer lock harvesting cannula
B 8FIC-V-169-LL	Type V 16ga X 9cm fat infiltration cannula
B 8FIC-I-187-LL	Type I 18ga X 7cm fat infiltration cannula
B 8FIC-I-195-LL	Type I 19ga X 5cm fat infiltration cannula
B 8FIC-I-203-LL	Type I 20ga X 3cm fat infiltration cannula
B 89LL3	Syringe-to-syringe luer lock-to-luer lock transfer, 3mm
B 89SR-10	Autoclavable syringe rack, 1cc - 10cc

LIPOSUCTION CANNULA SET

B 8MER-326L	Mercedes 3mm X 26cm large handle cannula
B 8MER-432L	Mercedes 4mm X 32cm large handle cannula
B 8BAS-432L	Basket 4mm X 32cm large handle cannula
B 8III-326L	Accelerator 3mm X 26cm large handle cannula
B 8III-432L	Accelerator 4mm X 32cm large handle cannula
B 8SPA-315L	Spatula 3mm X 15cm large handle cannula

TUMESCENT INFILTRATION SET

B 89055	Variable speed infiltration pump with volume read out
B 89IH-NC	Infiltration handle without controls
B 8TI-1232-LL	Tumescent infiltration cannula 12ga X 32cm
B 8TI-1426-LL	Tumescent infiltration cannula 14ga X 26cm
B 8TI-1415-LL	Tumescent infiltration cannula 14ga X 15cm
B 89213	Infiltration tubing for pump, 13', 10'/box

POWER LIPOSUCTION STARTER SET

B 8B-MER-326 PWR	Mercedes Power Cannula 3mm x 26cm
B 8B-MER-432 PWR	Mercedes Power Cannula 4mm x 32cm
B 8B-BAS-432 PWR	Basket Power Cannula 4mm x 32cm
B 8B-III-326 PWR	Accelerator Power Cannula 3mm x 26cm
B 8B-III-432 PWR	Accelerator Power Cannula 4mm x 32cm

VITRUVIAN™
Perfection in Body Sculpting

©2026 Black & Black Surgical

25-0525 2026 Ver. 02

5175 South Royal Atlanta Drive
Tucker, GA 30084
770.414.4880 phone
1.877.252.2517 toll free
770.414.4879 fax
info@blackandblacksurgical.com
www.blackandblacksurgical.com

